





# Micropipette Single Channel Autoclavable




Professional and Relieved

### PRODUCT DETAIL



### Identity Code


Press the ejection button to see the exclusive identity code of each pipette



Quality Guarantee

### Clear brand and specification

Joanlab or customized brand available  
Volume specification shows on each pipette help to identify



High Accuracy

### 3 Digits Volume Display


Clear readout and easy volume setting



Safe and Reliable

### Smooth operation

With ergonomic handle hook design for easy grip in the palm



Safe and Reliable

### CE Mark

The product has passed the CE certification of European Union,  
it can withstand high temperature sterilization of 121 Centigrade



## ACCESSORIES



① Calibration oil

(can not sent by express company as they refused)



② Tips



③ Accessory supplied



④ Quality control certificate



⑤ User manual

## PACKAGING




## MAIN FEATURES

**Light weight, ergonomic design.**  
**Digital display clearly reads volume setting.**  
**Wide measurement range, 0.1ul to 10000ul.**  
**With new material, it can be fully autoclavable.**  
**Ceramic piston, high corrosion resistance.**  
**PTFE sealing ring, more convenient to use.**

P Series Micropipette Autoclavable Technical Data						
Model	Volume Range	Testing Volume (µl)	Inaccuracy	Imprecision	Increment (µl)	Tip volume
P10	0.5 - 10ul	10	±1.00%	±1.00%	0.1ul	20ul
		5	±1.50%	±1.50%		
		1	±2.50%	±2.50%		
P20	2 - 20ul	20	±0.90%	±0.40%	0.5ul	200ul
		10	±1.20%	±1.00%		
P50	5 - 50ul	50	±0.60%	±0.30%	0.5ul	200ul
		25	±0.90%	±0.60%		
		5	±2.00%	±2.00%		
P100	10 - 100ul	100	±0.60%	±0.15%	1ul	200ul
		50	±0.90%	±0.40%		
		10	±2.00%	±0.10%		
P200	20 - 200ul	200	±0.60%	±0.15%	1ul	200ul
		100	±0.80%	±0.30%		
		20	±3.00%	±1.00%		
P1000	100 - 1000ul	1000	±0.60%	±0.20%	5ul	1000ul
		500	±0.70%	±0.25%		
		100	±1.50%	±0.70%		
P5000	1000 - 5000ul	5000	±0.50%	±0.15%	50ul	5000ul
		2500	±0.60%	±0.30%		
		1000	±0.70%	±0.30%		
P10000	2000 - 10000ul	10ml	±0.60%	±0.20%	0.1ml	10ml
		5ml	±1.20%	±0.30%		
		2ml	±1.50%	±0.40%		


**Authorized Distributor**  
**ALFA DIAGNOSTIK SDN BHD**  
**No 36, Taman Perindustrian Jaya, Jalan PJU 1A/13, Ara**  
**Damansara, 47301, Petaling Jaya, Selangor**  
**Phone: 03-78321000**  
**[www.alfadiagnostik.com.my](http://www.alfadiagnostik.com.my)**

# Micropipette Single Channel Half Autoclavable



Professional and Relieved

## PRODUCT DETAIL



## Identity Code

Each pipette has its own code, three-in-one for the code on the body, factory certificate and packaging box.




Quality Guarantee

High Accuracy

## Clear brand and specification

Joanlab or customized brand available  
Volume specification shows on each pipette help to identify



## Multiple specifications Wide measuring range

Capacity Range 0.1µl~10ml

Model: M2.5 \ M10 \ M20 \ M50 \ M100M200 \ M1000 \ M5000 \ M2000 \ M10000



Safe and Reliable

## Piano wire spring components

Good resilience, Resilience fast, no deformation



## ACCESSORIES



① Calibration oil

(can not sent by express company as they refused)



② Tips



③ Accessory supplied



④ Quality control certificate



⑤ User manual

## PACKAGING



## MAIN FEATURES

**Light weight, ergonomic design. Digital display clearly reads volume setting. Wide measurement rang, 0.1ul to 10000ul. With new material, it can be fully autoclavable. Ceramic piston, high corrosion resistance. PTFE sealing ring, more convenient to use.**

M Series Variable Volume Micropipette Technical Data						
Model	Volume Range	Testing Volume (µl)	Inaccuracy	Imprecision	Increment (µl)	Tips Volume
M2.5	0.1-2.5ul	1.25	±3.00%	±0.30%	0.05	10
M10	0.5-10ul	5	±1.50%	±1.50%	0.1	10
M20	2-20ul	10	±1.20%	±1.20%	0.5	200
M50	5-50ul	25	±0.90%	±0.30%	0.5	200
M100	10-100ul	50	±1.00%	±0.15%	1	200
M200	20-200ul	100	±0.70%	±0.15%	1	200
M1000	100-1000ul	500	±0.80%	±0.20%	5	1000
M5000	1000-5000ul	2500	±0.60%	±0.15%	50	5000
M10000	2000-10000ul	5ml	±1.20%	±0.20%	0.1ml	10ml

**Authorized Distributor**  
**ALFA DIAGNOSTIK SDN BHD**  
 No 36, Taman Perindustrian Jaya, Jalan PJU 1A/13, Ara  
 Damansara, 47301, Petaling Jaya, Selangor  
 Phone: 03-78321000  
[www.alfadiagnostik.com.my](http://www.alfadiagnostik.com.my)