

HEATING MANTLE



Heating Mantle



Intelligent temperature control



Heating evenly, heating up fast

TECHNICAL DATA

Model	HM		HMC		HMS			HMSC	
Capacity	100ml	250ml	500ml	1L	2L	3L	5L	10L	20L
Function	Heating		Heating		Heating/magnetic stirrer			Heating/magnetic stirrer	
Heating power	100W	150W	250W	350W	450W	600W	800W	1200W	2400W
Temperature control	/		Digital temperature control			/		Digital temperature control	
Max temperature	Room temp - 380°C		Room temp - 300°C		Room temp - 380°C			Room temp - 300°C	