

# INFRARED STERILIZER



## Infrared Sterilizer



Sterilization only takes  
5-7 seconds



Can be used in an anaerobic  
environment



No open-flame, more safe



Easy to clean

### TECHNICAL DATA

Model	DS-900	DS-900S
Power	150W	240W
Max. temperature inside the heating hole	825°C±50°C	
Heating temperature	Two temperatures available 400°C/850°C	
Max. diameter of sterilizer goods	φ14mm	φ35mm
Length of heating zone	150mm	100mm
Warm-up time	10 min	