

IMPROVACUTER® URINE COLLECTION SYSTEM

Streamline urine collection to improve diagnostic accuracy. From collection to analysis, IMPROVACUTER® Urine Products help you maximize the utility of one of the most collected and analyzed body fluids.

IMPROVACUTER® Urine Collection devices and evacuated tubes maintain a closed system designed to facilitate a simple and consistent, high-quality collection.

SCAN FOR
INSTRUCTIONS VIDEO



IMPROVACUTER® URINE CONTAINER CUP

IMPROVACUTER® Urine Collection Container is pre-attached a non-patient needle and provides a safer and easier alternative for specimen transfer. By using this closed system, pouring off, leakage and wet areas are avoided, which positively impacts the work environment and conditions.

- (120mL, PP material, yellow lid)

IMPROVACUTER® URINALYSIS TUBE

IMPROVACUTER® Urinalysis Tubes are used for urinary sediment and component analysis. The conical bottom is suitable for routine urine tests and microscopic sediment analysis. The preservative tube contain Sodium Propionate, Chlorhexidine and Ethyl p-hydroxybenzoate as additive, providing better sample preservation. Both with and without preservatives are available.

- without additive, 8mL/16x100mm, PET (conical base), conventional cap/yellow, paper label
- with preservatives, 8mL/16x100mm, PET (conical base), conventional cap/red, paper label

